



**Supplementary Figure 4. The 600K complex recognizes both Ub-GST and GST-Ub, but fails to ubiquitinate GST-Ub.** (A) The 600K complex ubiquitinates Ub-GST but not GST-Ub. Ubiquitination of Ub-GST (lane 1) and GST-Ub (lane 2) were examined as described in Figure 3A. (B) The 600K complex binds to Ub-GST and GST-Ub. The baculovirus-purified 600K complex was mixed with Ub-GST (lanes 1 and 4), GST-Ub (lanes 2 and 5) or GST (lanes 3 and 6) followed by pull down with Glutathione Sepharose. Pull down proteins (lanes 1-3) and 5% input (lanes 4-6) were separated by SDS-PAGE and subjected to immunoblotting with anti-HOIP, anti-HOIL-1L and anti-GST antibodies.